TECHNICAL BULLETIN

□NEW TEST □ TEST CHANGE

NOTIFICATION DATE: 06/05/2025 EFFECTIVE DATE: Immediate

Hepatitis C RNA Quant (HEPCRNA)

On the effective date, Northwell Health Laboratories will transition the Hepatitis C RT- PCR quantitative assay (HEPCRNA) from the Abbott Alinity m system to the Roche cobas 8800 system. The HCV test results and interpretations are comparable; the collection requirements are the same.

Test Requirements/Parameters	New	Current
Instrument	cobas 8800 System	Alinity m System
Methodology	RT- PCR	RT- PCR
Specimen Collection Container	Lavender-top (EDTA), Serum Separation Tubes (SST) or Plasma Preparation Tubes (PPT)	Lavender-top (EDTA), Serum Separator Tube (SST) or Plasma Preparation Tube (PPT)
Specimen Types	Plasma or Serum	Plasma or Serum
Specimen Volume /Minimum Volume	2.5 mL/0.65 mL (Note: Minimum volume does not allow for repeat testing) up to 6 days at 2-8°C	2.5 mL/0.75 mL (Note: Minimum volume does not allow for repeat testing) Up to 4 hours at 15-30°C
Specimen Storage	up to 12 weeks at \leq -18°C, up to 6 months at \leq -60°C	up to 3 days at 2-8°C up to 60 days at -20°C, Long term at ≤ -70°C
Specimen Transport	Refrigerated (2°C to 8°C)	Refrigerated (2°C to 8°C)
Reportable Range	15 to 100,000,000 IU/mL or 1.18 to 8 log ₁₀ IU/mL	12 to 100,000,000 IU/mL or 1.08 to 8 log ₁₀ IU/mL
Reference Range	Not Detected	Not Detected
Computer Interface Code	PDM# 6350030	PDM# 6350030
Test Order	HEPCRNA	HEPCRNA

If you have any questions, please contact Client Services at (800) 472-5757

