



Receiving AP Specimens in Cerner Millennium

POLICY: Anatomic pathology samples collected at Northwell Sites are ordered in the electronic medical record (EMR), Epic. Epic orders generate laboratory orders in Cerner Millennium and transferred to the laboratory. Receiving the specimens from the courier is reconciled electronically.

NOTE: Anatomic pathology samples collected at certain outreach locations may be ordered and requested on paper requisitions only.

PROCEDURE: This procedure describes the steps for receiving pathology specimens with electronic orders generated in Epic.

Step 1: In Cerner Millennium, click the Specimen Log-In By Accession icon on the app bar to open. Also open Pending Inquiry.	
Step 2: In Log-In by Accession. Accession should be chosen, then click Retrieve	
Step 3. Scan in each container. Choose Location at the bottom. Click Log-In	

Receiving AP Specimens in Cerner Millennium

Status changed from Dispatch to Collected and the Time Stamp for Received is Logged into Cerner

Log in By Accession

Task View Help

Pearl, Mary
DOB: 07/05/1950 NHPP Concorde... NHPP Concorde...
000547333
Age: 75 years Female

Collections Accession: 0 minutes ago

C	P	Accession	Coll Date	Coll Time	Coll ID	Order	Priority	Status	Coll Method
☐		99-24-163-15767	06/11/2024	14:11	AKOZER	Tissue Biopsy	RO-Routine Outreach	Collected	AP Specimen
☐		99-24-163-15768	06/11/2024	14:11	AKOZER	Tissue Biopsy	RO-Routine Outreach	Collected	AP Specimen

Location: 1111 Marcus - APS Log In Miss... Details

Ready P819 AKOZER 10/15/2025 14:57

Step 4. Continue on to
**Accessioning of Electronic
Anatomic Pathology
Specimens in Cerner
Millennium**