




Viewing Pending Pathology Orders in Cerner Millenium

POLICY: Anatomic pathology samples collected at Northwell Sites are ordered in the electronic medical record (EMR), Epic. Epic orders generate laboratory orders in Cerner Millenium. These orders can be viewed on a pending list prior to specimen receipt in the laboratory for specimen tracking and management.

PROCEDURE: Utilize the Cerner Millenium **Pending Inquiry** application to view pathology orders generated in Epic.

Step 1. Click the Pending Inquiry icon on the Cerner app bar.	
Step 2: a. The “Select Pending Procedures” box appears. b. Ensure <i>All Pending</i> is selected under STATUS. c. Click the magnifying glass under Test site to search for a collection location.	
Step 3. The “Test Site Lookup” window will appear. Find a location by typing in the initials of a collection location. For example, type PB for Peconic Bay and select <i>PB AP Grossing Bay</i> and click OK. - Pathology orders: “Prefix Grossing Bay” - Cytology orders: “Prefix CY Request”	
Step 3: This query generates a list of pending pathology orders that were collected in Epic but have not yet been received in the laboratory.	