



TECHNICAL BULLETIN

☐NEW TEST ☒TEST CHANGE

NOTIFICATION DATE: 06/09/2022

EFFECTIVE DATE: 06/20/2022

BK Virus DNA Test Update (BK Virus)

Northwell Health Laboratories would like to inform you that the BK Virus DNA (Plasma) test performed by Viracor-Eurofins will update the assay range from copies per mL to IU per mL.

Test Requirement/Parameters	New	Previous
Assay Range	33 to 3.30E+08 IU/mL Units: IU/mL (one IU is equal to 0.33 copies)	39 to 1.00E+10 copies/mL Units: copies/mL

Computer Interface Code: **PDM #5902135**

Test Order: **BK Virus**

If you have any questions, please contact Client Services at (516) 719-1100.