



## ***TECHNICAL BULLETIN***

**NEW TEST**    **TEST CHANGE**

NOTIFICATION DATE: MAY 14, 2021

EFFECTIVE DATE: May 18, 2021

### **Clozapine, Immunoassay (CLZ)**

On the effective date Northwell Health Laboratories will now offer Clozapine testing by Immunoassay as part of our continued quality improvement.

<b>Test Requirement/Parameters</b>	
<b>Method:</b>	Nanoparticle Homogenous Immunoassay
<b>Specimen Requirements:</b>	Red Top Tube with no additives or separator gel.
<b>Specimen Volume/Minimum Volume:</b>	500uL/200uL
<b>Specimen Stability:</b>	7 days refrigerated
<b>Reference Range:</b>	350-600 ng/mL
<b>Turnaround time</b>	24 hours

**CPT Coding: [80159](#)**

**Computer Interface Code: PDM [#2153767](#)**

**Test Order: [CLZ](#)**

**[If you have any questions, please contact Client Services at \(516\) 719-1100.](#)**