



TECHNICAL BULLETIN

☐ NEW TEST ☒ TEST CHANGE

NOTIFICATION DATE: 12/19/2019

EFFECTIVE DATE: 12/19/2019

Cytomegalovirus (CMV) by Qualitative PCR, Saliva CMVSAL

Explanation/Reason: On the effective date the *Cytomegalovirus (CMV) by qualitative PCR, Saliva test* will now be performed by Viracor Eurofins Clinical Diagnostics. The collection requirements will change from the special collection device ORA Collect OC-100 kit to Universal Transport Media (UTM). This change is a result of increasing difficulty in obtaining the special collection devices due to backorder from the current performing laboratory ARUP labs.

Test Requirement/Parameters	New Requirement	Previous Requirement
Specimen Requirements:	Universal Transport Media UTM	ORACollect OC-100 kit
Submission Container/Tube:	Universal Transport Media UTM	ORACollect OC-100 kit
Specimen Stability:	Room Temperature: 48 Hours Refrigerated or Frozen: 7 Days	Room Temperature or Refrigerated: 7 Days Frozen: 30 Days
Reference Range:	Not Detected	Negative

Computer Interface Code: **PDM # 1659370**

Test Order: **CMVSAL**

If you have any questions, please contact Client Services at (516) 719-1100.