



TECHNICAL BULLETIN

TEST UPDATE

NOTIFICATION DATE: 11/06/2020

EFFECTIVE DATE: 11/16/2020

Mycoplasma Pneumoniae IgM Antibody Reference Range and Interpretation Change (MYCO IGM)

Effective, November 16, 2020, the Reference Range and Interpretation for the Mycoplasma Pneumoniae IgM Antibody test will change as a result of the implementation of a new test reagent.

Test Parameters	New	Previous
Reference Range:	0.00-0.90 Index	<= 910 Units/mL
Interpretation:	Negative: <=0.90 Index Equivocal: 0.91-1.09 Index Positive: >=1.10 Index	Negative: <=910 Units/mL Low Positive: 911-1119 Units/mL Positive: >=1120 Units/mL

Computer Interface Code: [PDM 5700420](#)

Test Order: [MYCO IGM](#)

If you have any questions, please contact Client Services at (516) 719-1100.