



TECHNICAL BULLETIN

☐ NEW TEST ☒ TEST CHANGE

NOTIFICATION DATE: 08/16/2022

EFFECTIVE DATE: 08/22/2022

Tryptase (TRYPT)

On the effective date, Northwell Health Laboratories will now perform the Tryptase test on the Phadia 1000 instrument with improved turnaround time. This test will no longer be performed at LabCorp.

Test Requirements/Parameters	New	Previous
Methodology	ImmunoCAP Fluoroenzymeimmunoassay	ImmunoCAP Fluoroenzymeimmunoassay
Submission Container/Tube	Gold Top or Lavender Top Tube	Gold Top or Lavender Top Tube
Specimen Volume	0.7 mL	0.7 mL
Storage Instruction	Refrigerated	Refrigerated
Reference Range	<8.4 ug/L	2.2 - 13.2 ug/L
Computer Interface Code	PDM # 5913930	PDM # 5913930
Test Code	TRYPT	TRYPT

If you have any questions, please contact Client Services at (516) 719-1100.