

## TECHNICAL BULLETIN



## **⊠NEW TEST** □**TEST CHANGE**

NOTIFICATION DATE: 11/03/2022 EFFECTIVE DATE: Immediate

## Allergen-Specific IgE Testing - Stinging Insect Venom

Northwell Health Laboratories is pleased to announce the availability of component-resolved allergen-specific IgE tests for the following stinging insects: Honeybee (Apis mellifera), Paper Wasp (Polistes dominula), and Yellow Jacket/Common Wasp (Vespula vulgaris).

The advantages of component-resolved IgE testing for stinging insects include:

- Improved discrimination between species specific sensitization and cross-reactivity
- Improved identification of culprit venom(s)
- Increased likelihood of successful venom immunotherapy

<b>Test Parameters</b>	Yellow Jacket	Paper Wasp	<b>Honey Bee</b>	Venom
	<b>Component Panel</b>	<b>Component Panel</b>	<b>Component Panel</b>	<b>Component Panel</b>
Profile	rVes v 1 (i211), rVes	rPol d 5 (i210), CCD	rApi m 1 (i208), rApi	Honey Bee (i1), rApi m
Components	v 5 (i209), CCD	(o214), Tryptase	m 2 (i214), rApi m 3	1 (i208), rApi m 2
	(o214), Tryptase		(i215), rApi m 5	(i214), rApi m 3 (i215),
			(i216), rApi m 10	rApi m 5 (i216), rApi
			(i217),	m 10 (i217), rVes v 1
			CCD (o214),	(i211), rVes v 5 (i209),
			Tryptase	Paper Wasp (i4), rPol d
				5 (i210), Yellow hornet
				(i5), CCD, Tryptase
Method	ImmunoCAP	ImmunoCAP	ImmunoCAP	ImmunoCAP
	Fluorescence Enzyme	Fluorescence Enzyme	Fluorescence Enzyme	Fluorescence Enzyme
	Immunoassay	Immunoassay	Immunoassay	Immunoassay
Submission	Gold Top Tube	Gold Top Tube	Gold Top Tube	Gold Top Tube
Container/Tube				
Specimen	0.5 mL Serum	0.5 mL Serum	0.5 mL Serum	1.0 mL Serum
Volume				
Storage	Refrigerated	Refrigerated	Refrigerated	Refrigerated
Instruction				
Computer	PDM # 221782	PDM # 221784	PDM# 221783	PDM# 221781
Interface Code				
Test Code	RYELJACKCP	RPAPERWPCP	RHONEYBCP	VENOMCP
CPT Code	86008 x3, 83520	86008 x2, 83520	86008 x6, 83520	86008 x9, 86003 x3, 83520