

TECHNICAL BULLETIN

☐ NEW TEST ☒ TEST CHANGE

NOTIFICATION DATE: 11/22/2022

EFFECTIVE DATE: 11/30/2022

Drug Screen, Urine Routine Test Update (DRUGSCRUR)

On the Effective date, the Drug Screen, Urine Routine test will be updated. The Propoxyphene and the Methaqualone will be replaced with Oxycodone/Oxymorphone and Fentanyl.

Test Requirement/Parameters	New	Previous
Method	Screening test by Immunoassay (IA); confirmation of positive results by Mass Spectrometry (MS)	Gas Chromatography (GC) Liquid Chromatography/tandem mass spectrometry (LC/MS-MS)
Confirmatory Test Analytes	Oxycodone/Oxymorphone, Fentanyl, Amphetamines, Barbiturates, Benzodiazepines, Cannabinoid, Cocaine, Opiates, Phencyclidine, Methadone	Propoxyphene, Methaqualone, Amphetamines, Barbiturates, Benzodiazepines, Cannabinoid, Cocaine, Opiates, Phencyclidine, Methadone
Specimen Requirements	30 mL random urine. No preservative. Minimum: 3 – 5 mL	10-20 mL random urine. No preservative. Minimum: 5 – 10 mL
Specimen Stability	Room Temperature - 3 Days 14 Days Refrigerated	Room Temperature - 5 days

Computer Interface Code: **PDM #5940000**

Test Order: **DRUGSCRUR**

If you have any questions, please contact Client Services at (516) 719-1100.