

TECHNICAL BULLETIN

☐ NEW TEST ☒ TEST CHANGE

NOTIFICATION DATE: 04/30/2025

EFFECTIVE DATE: Immediately

Diluted Thrombin Time (DTT)

On the effective date, Northwell Health Laboratories will update the reference range for the Diluted Thrombin Time (DTT) assay to reflect a change in reagent dilution. This update improves the test's clinical utility, particularly for monitoring antithrombin therapy.

Test Requirements/Parameters	New	Current
Reference Range:	11 – 16.4 seconds	16 – 26 seconds
Computer Interface Code:	PDM# 5500555	PDM# 5500555
Test Order:	DTT	DTT



If you have any questions, please contact Client Services at (800) 472-5757.