TECHNICAL BULLETIN

⊠NEW TEST □**TEST CHANGE**

NOTIFICATION DATE: 10/05/2023 EFFECTIVE DATE: 10/10/2023

IgG Subsets (IGGSUB)

On the effective date, Northwell Health Laboratories, as part of its ongoing commitment to providing timely and accurate laboratory results, will now perform the **IgG Subsets** test, improving turnaround time. This test will no longer be performed at LabCorp.

IgG is the most abundant serum immunoglobulin in healthy individuals and is comprised of 4 subclasses (i.e., IgG1, IgG2, IgG3 and IgG4). Measurement of IgG subclasses is useful when evaluating patients with suspected humoral immunodeficiency.

Test Requirement	New	Previous
Specimen Type/Volume	3 mL serum (minimum 1.8 mL)	3 mL serum (minimum 1.8 mL)
Submission Container/Tube	Gold Top Tube (Serum separator tube)	Gold Top Tube (Serum separator tube)
Methodology	Quantitative Immunoturbidimetry	Quantitative Immunoturbidimetry
Transport Temperature	Refrigerated	Refrigerated
Computer Interface Code	PDM #5904190	PDM #5904190
Test Order	IGGSUB	IGGSUB