

☒ **NEW TEST** ☐ **TEST CHANGE**

NOTIFICATION DATE: 10/02/2023

EFFECTIVE DATE: 10/03/2023

Sesame Allergen Component (RSESAMECP)

On the effective date, Northwell Health Laboratories will now perform the Sesame Allergen Component test. This test supports the diagnosis of sesame allergy by detection of specific IgE antibodies to Ses i 1 (f449), which is one of the major allergens found in the sesame seeds (*Sesamum indicum*). Ses i 1 is reported to be the most clinically relevant allergen of sesame associated with severe allergic reactions, especially anaphylaxis.

Test Requirement	Sesame Allergen Component
Specimen Type/ Volume	0.8 mL serum (minimum 0.5 mL)
Submission Container/Tube	Gold top tube (Serum separator tube)
Methodology	ImmunoCAP Fluorescence Enzyme Immunoassay
Transport Temperature	Refrigerated
CPT Code	86008
Computer Interface Code	PDM #231707
Test Order	RSESAMECP



If you have any questions, please contact Client Services at (800) 472-5757.